



# Anti-STARD10 (aa 72-172) polyclonal antibody (DPABH-20195)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Antigen Description	May play specific roles in sperm maturation or fertilization.
Immunogen	antigen sequence, corresponding to amino acids 72-172 of Human PCTP-L. (Q9Y365)
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse, Rat, Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P, WB
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

## GENE INFORMATION

Gene Name	<a href="#">STARD10 StAR-related lipid transfer (START) domain containing 11 [ Homo sapiens ]</a>
-----------	---

<b>Official Symbol</b>	STARD10
<b>Synonyms</b>	STARD10; StAR-related lipid transfer (START) domain containing 10; PCTP2; CGI-52; NY-CO-28; SDCCAG28; PCTP-like protein; PCTP-L; antigen NY-CO-28; START domain containing 10; START domain-containing protein 10; stAR-related lipid transfer protein 10; serologically defined colon cancer antigen 28;
<b>Entrez Gene ID</b>	<a href="#">10809</a>
<b>Protein Refseq</b>	<a href="#">NP_006636.2</a>
<b>UniProt ID</b>	<a href="#">A0A024R5L8</a>
<b>Function</b>	lipid binding; protein binding;