



Anti-C1ORF64 (aa 91-167) polyclonal antibody (DPABH-20197)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	C1orf64 (chromosome 1 open reading frame 64) is a protein-coding gene. Diseases associated with C1orf64 include breast cancer.
Immunogen	antigen sequence: HLTRVTPMGG GCLAQARATL PLCRGSVASA SFPVSPLCPQ EVPEAKGKPV KAAPVRSSTW GTVKDSLKAL SSCVCGQ, corresponding to amino acids 91-167 of Human C1orf64. (Q8NEQ6)
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, IHC-P, WB
Format	Liquid
Size	100 µl
Buffer	pH: 7.20; Constituents: 59% PBS, 40% Glycerol
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	C1ORF64 chromosome 1 open reading frame 65 [Homo sapiens]
Official Symbol	C1ORF64
Synonyms	C1ORF64; chromosome 1 open reading frame 64; ERRF; uncharacterized protein C1orf64; ER-related factor;
Entrez Gene ID	149563
Protein Refseq	NP_849162.1
UniProt ID	Q8NEQ6