



Anti-BICC1 (aa 229-240) polyclonal antibody (DPABH-20484)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Putative RNA-binding protein. May be involved in regulating gene expression during embryonic development.
Immunogen	Synthetic peptide: C-QHISQTYNISVS, corresponding to internal sequence amino acids 229-240 of Human BICC1 (NP_001073981.1).
Isotype	IgG
Source/Host	Goat
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Format	Liquid
Size	100 µg
Buffer	Constituents: 0.5% BSA, Tris buffered saline, pH 7.3
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name [BICC1 BicC family RNA binding protein 2 \[Homo sapiens \]](#)

Official Symbol	BICC1
Synonyms	BICC1; BicC family RNA binding protein 1; BICC; CYSRD; protein bicaudal C homolog 1; bic-C; bicaudal C homolog 1; FGFR2-BICC1 fusion kinase protein;
Entrez Gene ID	80114
Protein Refseq	NP_001073981.1
UniProt ID	Q9H694
Function	poly(A) RNA binding;