



Anti-TMEM30A (aa 108-157) polyclonal antibody (DPABH-21088)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Belongs to the CDC50/LEM3 family.
Immunogen	Synthetic peptide corresponding to a region within N terminal amino acids 108-157 (FTLEKSFEEN VFMYYGLSNF YQNHRRYVKS RDDSQLNGDS SALLNPSKEC) of Human Transmembrane protein 30A (NP_060717).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, IHC-P
Format	Liquid
Size	50 µg
Buffer	Constituents: 2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name	TMEM30A transmembrane protein 31A [Homo sapiens]
Official Symbol	TMEM30A
Synonyms	TMEM30A; transmembrane protein 30A; CDC50A; C6orf67; cell cycle control protein 50A;
Entrez Gene ID	55754
Protein Refseq	NP_001137430.1
UniProt ID	Q9NV96
Function	protein binding;