



Anti-TAZ (aa 100-200) polyclonal antibody (DPABH-21226)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Some isoforms may be involved in cardiolipin (CL) metabolism.
Immunogen	Synthetic peptide conjugated to KLH derived from within residues 100 - 200 of Human Tafazzin/ TAZ.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse, Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Size	100 µg
Buffer	pH: 7.40; Constituent: PBS
Preservative	0.02% Sodium Azide
Storage	Store at 4°C short term (1-2 weeks). Aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name	TAZ tafazzin [Homo sapiens]
Official Symbol	TAZ
Synonyms	TAZ; tafazzin; EFE; BTHS; EFE2; G4.5; Taz1; CMD3A; LVNCX; protein G4.5;
Entrez Gene ID	6901
Protein Refseq	NP_000107.1
UniProt ID	Q16635
Pathway	Acyl chain remodeling of CL; Glycerophospholipid metabolism; Metabolism; Metabolism of proteins.
Function	1-acylglycerophosphocholine O-acyltransferase activity;