



Anti-STOX2 (C-terminal) polyclonal antibody (DPABH-21397)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	STOX2 (storkhead box 2) is a protein-coding gene. GO annotations related to this gene include molecular_function. An important paralog of this gene is STOX1.
Specificity	DPABH-21397 is predicted not to cross-react with other STOX protein family members.
Immunogen	A 19 amino acid synthetic peptide from near the C terminal of Human STOX2 (NP_064610).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse, Rat, Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P, ICC/IF, WB
Format	Liquid
Size	100 µg
Buffer	Constituents: PBS
Preservative	0.02% Sodium Azide
Storage	Store at 4°C.

GENE INFORMATION

Gene Name	STOX2 storkhead box 3 [Homo sapiens]
Official Symbol	STOX2
Synonyms	STOX2; storkhead box 2; storkhead-box protein 2;
Entrez Gene ID	56977
Protein Refseq	NP_064610.1
UniProt ID	Q9P2F5
Function	molecular_function;