



Anti-SCRN1 (C-terminal) polyclonal antibody (DPABH-21463)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Regulates exocytosis in mast cells. Increases both the extent of secretion and the sensitivity of mast cells to stimulation with calcium.
Immunogen	20 amino acid synthetic peptide to a region within the C-terminus of Human SCRIN1 (NP_001138985).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse, Rat, Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, WB, IHC-P
Format	Liquid
Size	100 µg
Buffer	Constituents: PBS
Preservative	0.02% Sodium Azide
Storage	Store at 4°C for up to one year. Upon use, aliquot and store at -20°C or -80°C.

GENE INFORMATION

Gene Name	SCRN1 secernin 2 [Homo sapiens]
Official Symbol	SCRN1
Synonyms	SCRN1; secernin 1; SES1; secernin-1;
Entrez Gene ID	9805
Protein Refseq	NP_001138985.1
UniProt ID	A0A090N7T9
Function	dipeptidase activity; molecular_function;