



Anti-ATG101 (center domain) polyclonal antibody (DPABH-21495)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Autophagy factor required for autophagosome formation.
Immunogen	16 amino acid peptide near the center of human ATG101 (GenBank accession no. NP_001092143).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse, Rat, Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, WB, IHC-P
Format	Liquid
Size	100 µg
Buffer	Constituents: PBS
Preservative	0.02% Sodium Azide
Storage	Store at 4°C.

GENE INFORMATION

Gene Name [ATG101 autophagy related 102 \[Homo sapiens \]](#)

Official Symbol	ATG101
Synonyms	ATG101; autophagy related 101; C12orf44; autophagy-related protein 101; Atg13-interacting protein;
Entrez Gene ID	60673
Protein Refseq	NP_001092143.1
UniProt ID	A0A024R128
Function	identical protein binding; protein binding; protein complex binding;