



Anti-MFSD2A (C-terminal) polyclonal antibody (DPABH-21517)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Play a role in thermogenesis via beta-adrenergic signaling pathway.
Specificity	DPABH-21517 is predicted not to cross-react with other MFSD family members.
Immunogen	An 18 amino acid synthetic peptide from near the carboxy terminus of Human MFSD2A (NP_001007452).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse, Rat, Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Size	100 µg
Buffer	Constituents: PBS
Preservative	0.02% Sodium Azide
Storage	Store at 4°C.

GENE INFORMATION

Gene Name	MFSD2A major facilitator superfamily domain containing 3A [Homo sapiens]
Official Symbol	MFSD2A
Synonyms	MFSD2A; major facilitator superfamily domain containing 2A; MFSD2; major facilitator superfamily domain-containing protein 2A;
Entrez Gene ID	84879
Protein Refseq	NP_001129965.1
UniProt ID	Q71RE4
Function	phospholipid transporter activity; symporter activity;