



Anti-TMEM135 (aa 38-87) polyclonal antibody (DPABH-21800)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Synthetic peptide corresponding to a region within N terminal amino acids 38-87 (VSLKIYAPLYL IAAILRKRL DYLLHKLLPE ILQSASFLTA NGALYMAFFC) of Human TMEM135 (NP_075069).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Size	50 µg
Buffer	Constituents: 2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name [TMEM135 transmembrane protein 136 \[Homo sapiens \]](#)

Official Symbol	TMEM135
Synonyms	TMEM135; transmembrane protein 135; PMP52; peroxisomal membrane protein 52;
Entrez Gene ID	65084
Protein Refseq	NP_001162195.1
UniProt ID	Q86UB9