



Anti-PORCN (aa 38-87) polyclonal antibody (DPABH-21827)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Modulates the processing of Wnt proteins. Probable protein-cysteine N-palmitoyltransferase that palmitoylates Wnt family members.
Immunogen	Synthetic peptide corresponding to a region within N terminal amino acids 38-87 (ACRLLWRLGL PSYKHASTV AGGFFSLYHF FQLHMVWVVL LSLLCYLVLF) of Human PORCN (NP_035050)
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, ICC/IF
Format	Liquid
Size	50 µg
Buffer	Constituents: 2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	PORCN porcupine homolog (Drosophila) [Homo sapiens]
Official Symbol	PORCN
Synonyms	PORCN; porcupine homolog (Drosophila); PPN; DHOF; FODH; MG61; PORC; protein-cysteine N-palmitoyltransferase porcupine; probable protein-cysteine N-palmitoyltransferase porcupine;
Entrez Gene ID	64840
Protein Refseq	NP_001269096.1
UniProt ID	Q9H237
Pathway	Signal Transduction; WNT ligand biogenesis and trafficking; Wnt signaling pathway;
Function	transferase activity, transferring acyl groups;