



Rabbit Anti-Human LPPR2 Polyclonal antibody (DPABH-21890)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	LPPR2 fusion protein, sequence: TGAAIATFLVTCVVHNFQSRPPSGRRLSPWEDLGQAPTMDSPLEKNPRSAGRIRHRHGSP HPSRRTAPAVAT (272-343aa encoded by BC009378)
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human, Mouse
Purification	Antigen affinity purification
Conjugate	Unconjugated
Applications	WB, IHC, ELISA
Positive Control	3T3-L1 cells, HeLa cells
Format	Liquid
Size	50 µl, 100 µl
Buffer	PBS with 0.1% sodium azide and 50% glycerol pH 7.3.
Preservative	0.1% Sodium Azide
Storage	Store at -20°C. Aliquoting is unnecessary for -20°C storage.

BACKGROUND

Keywords LPPR2; lipid phosphate phosphatase-related protein type 2; PRG4; PRG-4; plasticity related gene 4; plasticity-related gene 4 protein;

GENE INFORMATION

Entrez Gene ID [64748](#)

UniProt ID [Q96GM1](#)
