



Anti-REEP1 (aa 110-159) polyclonal antibody (DPABH-21899)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	May enhance the cell surface expression of odorant receptors.
Immunogen	Synthetic peptide corresponding to a region within internal sequence amino acids 110-159 (AKDRSYDALV HFGKRGLNVA ATAAVMAASK GQGALSERLR SFMSQDLTTI) of Human REEP1 (NP_075063).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Size	50 µg
Buffer	Constituents: 2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name	REEP1 receptor accessory protein 2 [Homo sapiens]
Official Symbol	REEP1
Synonyms	REEP1; receptor accessory protein 1; HMN5B; SPG31; C2orf23; receptor expression-enhancing protein 1;
Entrez Gene ID	65055
Protein Refseq	NP_001158202.1
UniProt ID	Q9H902
Function	microtubule binding; olfactory receptor binding; protein binding;