



Anti-TXNDC15 (aa 254-303) polyclonal antibody (DPABH-21970)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	TXNDC15 (thioredoxin domain containing 15) is a protein-coding gene.
Immunogen	Synthetic peptide corresponding to a region within C terminal amino acids 254-303 (STRFGTVAVP NILLFQGAKP MARFNHTDRT LETLKIFIFN QTGIEAKKNV) of Human TXNDC15 (NP_078991).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Size	50 µg
Buffer	Constituents: 2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name	TXNDC15 thioredoxin domain containing 16 [Homo sapiens]
Official Symbol	TXNDC15
Synonyms	TXNDC15; thioredoxin domain containing 15; UNQ335; C5orf14; thioredoxin domain-containing protein 15; 2310047H23Rik; disulfide isomerase;
Entrez Gene ID	79770
Protein Refseq	NP_078991.3
UniProt ID	Q7Z345