



Anti-FHOD3 (aa 390-439) polyclonal antibody (DPABH-21992)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Actin-organizing protein that may cause stress fiber formation together with cell elongation.
Immunogen	Synthetic peptide corresponding to a region within C terminal amino acids 390 - 439 (SGKFSGSSPA PPSQPQGLSY AEDAAEHENM KAVLKTSSPS VEDATPALGV) of Human FHOD3 (AAH81563).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse, Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, IHC-P
Format	Liquid
Size	50 µg
Buffer	Constituents: 2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name	FHOD3 formin homology 2 domain containing 4 [Homo sapiens]
Official Symbol	FHOD3
Synonyms	FHOD3; formin homology 2 domain containing 3; FHOS2; Formactin2; FH1/FH2 domain-containing protein 3; hFHOS2; formactin-2; formin homolog overexpressed in spleen 2;
Entrez Gene ID	80206
Protein Refseq	NP_001268668.1
UniProt ID	Q2V2M9
Function	actin binding; protein binding;