



Anti-SP6 (aa 327-376) polyclonal antibody (DPABH-22003)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Promotes cell proliferation.
Immunogen	Synthetic peptide corresponding to a region within C terminal amino acids 327 - 376 (DHLAKHMKTH EGAKEEAAGA ASGEGKAGGA VEPPGGKGKR EAEGSVAPSN) of Human SP6 (NP_954871).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Rat, Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P, WB
Format	Liquid
Size	50 µg
Buffer	Constituents: 2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name	SP6 Sp7 transcription factor [Homo sapiens]
Official Symbol	SP6
Synonyms	SP6; Sp6 transcription factor; EPFN; KLF14; EPIPROFIN; transcription factor Sp6; Kruppel-like factor 14; krueppel-like factor 14;
Entrez Gene ID	80320
Protein Refseq	NP_001245177.1
UniProt ID	B3KXP2
Function	DNA binding; metal ion binding;