



Anti-ZNF93 (aa 182-231) polyclonal antibody (DPABH-22012)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	May be involved in transcriptional regulation.
Immunogen	Synthetic peptide corresponding to a region within N terminal amino acids 182-231 (FNQFSTLITH KKIHTGEKPY ICEECGKAFK YSSALNTHKR IHTGEKPYKC) of Human ZNF93 (NP_112495).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Size	50 µg
Buffer	Constituents: 2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name	ZNF93 zinc finger protein 94 [Homo sapiens]
Official Symbol	ZNF93
Synonyms	ZNF93; zinc finger protein 93; TF34; HPF34; HTF34; ZNF505; zinc finger protein 505; zinc finger protein HTF34;
Entrez Gene ID	81931
Protein Refseq	NP_112495.2
UniProt ID	P35789
Function	DNA binding; sequence-specific DNA binding transcription factor activity; zinc ion binding;