



Anti-CRISPLD2 (aa 1-50) polyclonal antibody (DPABH-22034)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Promotes matrix assembly.
Immunogen	Synthetic peptide corresponding to a region within N terminal amino acids 1-50 (MSCVLGGVIP LGLLFLVCGS QGYLLPNVTL LEELLSKYQH NESH SRVRRRA) of Human CRISPLD2 (NP_113664).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Size	50 µg
Buffer	Constituents: 2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name	CRISPLD2 cysteine-rich secretory protein LCCL domain containing 3 [Homo sapiens]
Official Symbol	CRISPLD2
Synonyms	CRISPLD2; cysteine-rich secretory protein LCCL domain containing 2; CRISP11; LCRISP2; cysteine-rich secretory protein LCCL domain-containing 2; CRISP-11; cysteine-rich secretory protein 11; LCCL domain containing cysteine-rich secretory protein 2; LCCL domain-containing cysteine-rich secretory protein 2;
Entrez Gene ID	83716
Protein Refseq	NP_113664.1
UniProt ID	Q9H0B8
Function	heparin binding;