



Anti-SPNS1 (aa 2-51) polyclonal antibody (DPABH-22045)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Sphingolipid transporter (By similarity). May be involved in necrotic or autophagic cell death.
Immunogen	Synthetic peptide corresponding to a region within N terminal amino acids 2-51 (AGSDTAPFLS QADDPDDGPV PGTPGLPGST GNPKSEEPEV PDQEGLQRIT) of Human SPNS1 (NP_114427).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Size	50 µg
Buffer	Constituents: 2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name	SPNS1 spinster homolog 2 (Drosophila) [Homo sapiens]
Official Symbol	SPNS1
Synonyms	SPNS1; spinster homolog 1 (Drosophila); LAT; nrs; SPIN1; SPINL; HSpin1; PP2030; protein spinster homolog 1; spinster-like protein 1;
Entrez Gene ID	83985
Protein Refseq	NP_001135920.1
UniProt ID	Q9H2V7
Function	protein binding;