



Anti-ZDHHC16 (aa 111-160) polyclonal antibody (DPABH-22056)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	May be involved in apoptosis regulation.
Immunogen	Synthetic peptide corresponding to a region within N terminal amino acids 111-160 (SVPRLCWHFF YSHWNLILIV FHYYQAITTP PGYPPQGRND IATVSICKKC) of Human ZDHHC16 (NP_115703).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Size	50 µg
Buffer	Constituents: 2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name	ZDHHC16 zinc finger, DHHC-type containing 17 [Homo sapiens]
Official Symbol	ZDHHC16
Synonyms	ZDHHC16; zinc finger, DHHC-type containing 16; APH2; DHHC-16; probable palmitoyltransferase ZDHHC16; Abl-philin 2; zinc finger, DHHC domain containing 16; zinc finger DHHC domain-containing protein 16;
Entrez Gene ID	84287
Protein Refseq	NP_001274732.1
UniProt ID	B4DNL2
Function	palmitoyltransferase activity; protein-cysteine S-palmitoyltransferase activity; zinc ion binding;