



Anti-ZNF512 (aa 205-254) polyclonal antibody (DPABH-22066)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	May be involved in transcriptional regulation.
Immunogen	Synthetic peptide corresponding to a region within internal sequence amino acids 205-254 (QLRSLAGMKY HVMANHNLSLP ILKAGDEIDE PSERERLRTV LKRLGKLRCM) of Human ZNF512 (NP_115810).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Size	50 µg
Buffer	Constituents: 2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name	ZNF512 zinc finger protein 513 [Homo sapiens]
Official Symbol	ZNF512
Synonyms	ZNF512; zinc finger protein 512;
Entrez Gene ID	84450
Protein Refseq	NP_001258215.1
UniProt ID	Q96ME7
Function	DNA binding; metal ion binding;