



Anti-TMEM69 (aa 131-180) polyclonal antibody (DPABH-22102)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	TMEM69 (transmembrane protein 69) is a protein-coding gene.
Immunogen	Synthetic peptide corresponding to a region within internal amino acids 131-180 (AYGASFLSFLGGIRWGFALPEGSPA KPDYLN LASSAAPLFFSWFAFLIS E) of Human TMEM69 (NP_057570).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, IHC-P
Format	Liquid
Size	50 µg
Buffer	Constituents: 2% Sucrose, PBS
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	TMEM69 transmembrane protein 70 [Homo sapiens]
Official Symbol	TMEM69
Synonyms	TMEM69; transmembrane protein 69; C1orf154;
Entrez Gene ID	51249
Protein Refseq	NP_057570.2
UniProt ID	Q5SWH9