



Anti-CETN1 (aa 1-162) polyclonal antibody (DPABH-22878)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Plays a fundamental role in microtubule-organizing center structure and function.
Immunogen	Recombinant protein fragment containing a sequence corresponding to a region within amino acids 1-162 of Human Centrin 1 (AAH29515).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Size	50 µl
Buffer	Constituents: 20% Glycerol, 0.1M Tris, 0.1M Glycine, pH 7
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name	CETN1 centrin, EF-hand protein, 2 [Homo sapiens]
Official Symbol	CETN1
Synonyms	CETN1; centrin, EF-hand protein, 1; CEN1; CETN; centrin-1; caltractin; EF-hand protein; calcium binding protein;
Entrez Gene ID	1068
Protein Refseq	NP_004057.1
UniProt ID	Q12798
Function	ATP binding; ATP-dependent helicase activity; G-protein beta/gamma-subunit complex binding; calcium ion binding