



Rabbit Anti-Human Bcl2 Polyclonal Antibody (DPABH-22974)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	A synthetic peptide of human BCL2
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human, Mouse, Rat
Purification	Affinity purification
Conjugate	Unconjugated
Applications	WB, IHC, IF
Positive Control	Raji, HL-60, U-937
Format	Liquid
Size	50 µl, 100 µl
Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Preservative	0.02% Sodium Azide
Storage	Store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Introduction	This gene encodes an integral outer mitochondrial membrane protein that blocks the apoptotic death of some cells such as lymphocytes. Constitutive expression of BCL2, such as in the case
---------------------	--

of translocation of BCL2 to Ig heavy chain locus, is thought to be the cause of follicular lymphoma. Alternative splicing results in multiple transcript variants.

Keywords	BCL2; B-cell CLL/lymphoma 2; Bcl-2; PPP1R50; apoptosis regulator Bcl-2; protein phosphatase 1, regulatory subunit 50;
-----------------	---

GENE INFORMATION

Entrez Gene ID	596
-----------------------	---------------------

UniProt ID	P10415
-------------------	------------------------
