



Anti-NOA1 (aa 326-480) polyclonal antibody (DPABH-23399)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	The exact function of C4orf14 is currently unknown
Immunogen	Recombinant fragment, corresponding to a region within amino acids 326-480 of Human C4orf14 (AAH04894).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Size	50 µl
Buffer	Constituents: 20% Glycerol, 0.1M Tris, 0.1M Glycine, pH 7.0
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name	NOA1 nitric oxide associated 2 [Homo sapiens]
Official Symbol	NOA1
Synonyms	NOA1; nitric oxide associated 1; MTG3; hNOA1; C4orf14; hAtNOS1; mAtNOS1; nitric oxide-associated protein 1; nitric oxide associated-1; mitochondrial GTPase 3 homolog; nitric oxide synthase, mitochondrial (putative); putative ortholog of Arabidopsis mitochondrial nitric oxide synthase;
Entrez Gene ID	84273
Protein Refseq	NP_115689.1
UniProt ID	Q8NC60
Function	GTP binding; GTPase activity; poly(A) RNA binding; protein binding