



Anti-LRFN2 (C-terminal) polyclonal antibody (DPABH-23481)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Promotes neurite outgrowth in hippocampal neurons. Enhances the cell surface expression of 2 NMDA receptor subunits GRIN1 and GRIN2A. May play a role in redistributing DLG4 to the cell periphery.
Immunogen	Synthetic peptide corresponding to 14 aa from near the C-terminus of Human LRFN2 (Q9ULH4).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse, Rat, Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, IHC-P, WB
Format	Liquid
Size	100 µg
Buffer	Constituents: PBS
Preservative	0.02% Sodium Azide
Storage	Store at 4°C, stable for one year. Should not be exposed to prolonged high temperatures. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	LRFN2 leucine rich repeat and fibronectin type III domain containing 3 [Homo sapiens]
Official Symbol	LRFN2
Synonyms	LRFN2; leucine rich repeat and fibronectin type III domain containing 2; SALM1; FIGLER2; KIAA1246; leucine-rich repeat and fibronectin type-III domain-containing protein 2; synaptic adhesion-like molecule 1; fibronectin type III, immunoglobulin and leucine rich repeat domains 2;
Entrez Gene ID	57497
Protein Refseq	NP_065788.1
UniProt ID	Q9ULH4