



Anti-LRFN4 (C-terminal) polyclonal antibody (DPABH-23490)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Promotes neurite outgrowth in hippocampal neurons. May play a role in redistributing DLG4 to the cell periphery.
Immunogen	Synthetic peptide corresponding to 14 aa from near the C-terminus of Human LRFN4 (Q6PJG9).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse, Rat, Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC/IF, IHC-P, WB
Format	Liquid
Size	100 µg
Buffer	Constituents: PBS
Preservative	0.02% Sodium Azide
Storage	Store at 4°C, stable for one year. Should not be exposed to prolonged high temperatures. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	LRFN4 leucine rich repeat and fibronectin type III domain containing 5 [Homo sapiens]
Official Symbol	LRFN4
Synonyms	LRFN4; leucine rich repeat and fibronectin type III domain containing 4; SALM3; SALM3.; FIGLER6; leucine-rich repeat and fibronectin type-III domain-containing protein 4; fibronectin type III, immunoglobulin and leucine rich repeat domains 6;
Entrez Gene ID	78999
Protein Refseq	NP_076941.2
UniProt ID	A0A024R5I6