



Rabbit Anti-Human IRGC Polyclonal antibody (DPABH-23677)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	IRGC fusion protein, sequence: MATSKLPVVPGEEENTILMAKERLEALRTAFESGDLPPQAASHLQELLASTESIRLEVGV GESGAGKSSLINALRGLEAEDPGAALTGVMETTMQPSYPHPQFPDVTLWDLPGAGSPGC PADKYLKQVDFSRDYDFLLVSPRRCGAVETRLAAEILCQGKKFYFVRTKVDLDAATRTQ RPSGFREAAVLQEIRDHCAERLREAGVADPRIFLVSNLSPARYDFPTLVSTWEHDLP SHR RHAGLLSLPDISLEALQKKKAMLQEQLKTALVLGVIQALPVPGLAAAYDDALLI HSLRG YHRSFGLDDDSLAKLAEQVGKQAGDLRSVIRSPANEVSPETVLR LYSQSSDGAMRVARA FERGIPVFGTLVAGGISFGAVYTMLQGCLNEMAEDAQRVRI KALEDDEPQPEVSLEVASD NGVEKGGSGEGGGEEAPLSTCRKLGLLLKYILDSWKKH DSEEK (1-46)
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human, Mouse, Rat
Purification	Antigen affinity purification
Conjugate	Unconjugated
Applications	WB, ELISA
Positive Control	mouse testis tissue
Format	Liquid
Size	50 µl, 100 µl
Buffer	PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Preservative	0.02% Sodium Azide

Storage

Store at -20°C. Aliquoting is unnecessary for -20°C storage.

BACKGROUND

Keywords

IRGC; immunity-related GTPase family, cinema; IFGGE; IRGC1; ligp5; CINEMA; R30953_1; interferon-inducible GTPase 5; interferon inducible GTPase 5; immunity-related GTPase cinema 1; immunity-related GTPase family, cinema 1; interferon-gamma-inducible GTPase IFGGE protein;

GENE INFORMATION

Entrez Gene ID

[56269](#)

UniProt ID

[J7NNX4](#)
