



Anti-ARFGAP2 (aa 63-328) polyclonal antibody (DPABH-23865)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	GTPase-activating protein (GAP) for ADP ribosylation factor 1 (ARF1). Implicated in coatamer-mediated protein transport between the Golgi complex and the endoplasmic reticulum. Hydrolysis of ARF1-bound GTP may lead to dissociation of coatamer from Golgi-derived membranes to allow fusion with target membranes.
Immunogen	Recombinant protein fragment containing a sequence corresponding to a region within amino acids 63-328 of Human ZNF289.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Size	50 µl
Buffer	Constituents: 20% Glycerol, 1% BSA, 1X PBS, pH 7.0
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name	ARFGAP2 ADP-ribosylation factor GTPase activating protein 3 [Homo sapiens]
Official Symbol	ARFGAP2
Synonyms	ARFGAP2; ADP-ribosylation factor GTPase activating protein 2; IRZ; ZFP289; ZNF289; NBLA10535; ADP-ribosylation factor GTPase-activating protein 2; GTPase-activating protein ZNF289; zinc finger protein 289, ID1 regulated;
Entrez Gene ID	84364
Protein Refseq	NP_001229761.1
UniProt ID	G5E9L0
Pathway	Endocytosis;
Function	ARF GTPase activator activity; zinc ion binding;