



Anti-LZTR1 (N-terminal) polyclonal antibody (DPABH-24057)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Probable transcriptional regulator that may play a crucial role in embryogenesis.
Specificity	At least four isoforms of LZTR1 are known to exist; DPABH-24057 will recognize all but the shortest isoform. DPABH-24057 is predicted not to cross-react with LZTR2.
Immunogen	A 14 amino acid synthetic peptide near the N terminus of Human LZTR1 (CAJ86451).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse, Rat, Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P, WB, ICC/IF
Format	Liquid
Size	100 µg
Buffer	Constituents: PBS
Preservative	0.02% Sodium Azide
Storage	Store at 4°C, stable for one year. Should not be exposed to prolonged high temperatures. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	LZTR1 leucine-zipper-like transcription regulator 2 [Homo sapiens]
Official Symbol	LZTR1
Synonyms	LZTR1; leucine-zipper-like transcription regulator 1; BTBD29; LZTR-1; SWNTS2; leucine-zipper-like transcriptional regulator 1; leucine-zipper-like regulator-1;
Entrez Gene ID	8216
Protein Refseq	NP_006758.2
UniProt ID	Q8N653
Function	sequence-specific DNA binding transcription factor activity;