



Anti-LYRM2 (C-terminal) polyclonal antibody (DPABH-24132)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Specificity	DPABH-24132 is predicted not to crossreact with other LYRM family members.
Immunogen	A 16 amino acid peptide near the carboxy terminus of Human LYRM2 (NP_065199).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse, Rat, Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, ICC, ICC/IF
Format	Liquid
Size	100 µg
Buffer	Constituents: PBS
Preservative	0.02% Sodium Azide
Storage	Store at 4°C, stable for one year. Should not be exposed to prolonged high temperatures. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name [LYRM2 LYR motif containing 3 \[Homo sapiens \]](#)

Official Symbol	LYRM2
Synonyms	LYRM2; LYR motif containing 2; DJ122O8.2; LYR motif-containing protein 2;
Entrez Gene ID	57226
Protein Refseq	NP_065199.1
UniProt ID	Q9NU23