



Anti-TMEM184C (C-terminal) polyclonal antibody (DPABH-24188)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Possible tumor suppressor which may play a role in cell growth.
Specificity	DPABH-24188 predicted to not crossreact with other members of the TMEM184 protein family.
Immunogen	A 15 amino acid peptide near the carboxy terminus of Human TMEM184C (NP_060711).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse, Rat, Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	IHC-P, WB, ICC/IF
Format	Liquid
Size	100 µg
Buffer	Constituents: PBS
Preservative	0.02% Sodium Azide
Storage	Store at 4°C, stable for one year. Should not be exposed to prolonged high temperatures. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	TMEM184C transmembrane protein 185C [Homo sapiens]
Official Symbol	TMEM184C
Synonyms	TMEM184C; transmembrane protein 184C; TMEM34; transmembrane protein 34;
Entrez Gene ID	55751
Protein Refseq	NP_060711.2
UniProt ID	Q9NVA4