



Rabbit Anti-Human PHYHIPL Polyclonal antibody (DPABH-24298)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	PHYHIPL fusion protein, sequence: MEVPRLDHALNSPTSPCEEVIKNLSLEAIQLCDRDGNKSQDSGIAEMEELPVPHNIKISN ITCDSFKISWEMDSKSKDRITHYFIDLNNKENKNSNFKHKHDVPTKLVAKAVPLPMTVRG HWFLSPRTEYTVAVQTASKQVDGDYVSEWSEIEFCTADYSKVHLTQLLEKAEVIAGRM LKFSVFYRNQHKEYFDYVREHHGNAMQPSVKDNSGSHGSPISGKLEGIFFSCSTEFNTGK PPQDSPYGRYRFEIAAEKLFNPNTNLYFGDFYCMYTAYHYVILVIAPVGSPGDEFCKQRL PQLNSKDNKFLTCTEEDGVLVYHHAQDVILEVIYTDVPDLSLGTVAEITGHQLMSLSTAN AKKDPSCKTCNISVGR (1-376aa encoded by BC011268)
------------------	---

Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human, Mouse, Rat
Purification	Antigen affinity purification
Conjugate	Unconjugated
Applications	WB, ELISA
Positive Control	mouse kidney tissue, mouse brain tissue, mouse testis tissue
Format	Liquid
Size	50 µl, 100 µl
Buffer	PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Preservative	0.02% Sodium Azide

Storage

Store at -20°C. Aliquoting is unnecessary for -20°C storage.

BACKGROUND

Keywords

PHYHIPL; phytanoyl-CoA 2-hydroxylase interacting protein-like; phytanoyl-CoA hydroxylase-interacting protein-like;

GENE INFORMATION

Entrez Gene ID

[84457](#)

UniProt ID

[Q96FC7](#)
