



Anti-THRB (aa 201-212) polyclonal antibody (DPABH-24434)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	High affinity receptor for triiodothyronine.
Specificity	DPABH-24434 is not expected to cross-react to THRA.
Immunogen	Synthetic peptide: C-REELQKSIGHKP, corresponding to internal sequence amino acids 201-212 of Human Thyroid Hormone Receptor beta (NP_000452.2).
Isotype	IgG
Source/Host	Goat
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Size	100 µg
Buffer	Constituents: 0.5% BSA, Tris buffered saline, pH 7.3
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	THRB thyroid hormone receptor, beta [Homo sapiens]
Official Symbol	THRB
Synonyms	THRB; thyroid hormone receptor, beta; GRTH; PPTH; THR1; ERBA2; NR1A2; THRB1; THRB2; C-ERBA-2; C-ERBA-BETA; thyroid hormone receptor beta; oncogene ERBA2; nuclear receptor subfamily 1 group A member 2; thyroid hormone nuclear receptor beta variant 1; thyroid hormone receptor, beta (erythroblastic leukemia viral (v-erb-a) oncogene homolog 2, avian);
Entrez Gene ID	7068
Protein Refseq	NP_000452.2
UniProt ID	F1D8N7
Pathway	Gene Expression; Neuroactive ligand-receptor interaction; Nuclear Receptor transcription pathway; SIDS Susceptibility Pathways
Function	DNA binding; chromatin DNA binding; enzyme binding; protein binding