



Anti-ANKRD16 (aa 287-316) polyclonal antibody (DPABH-24832)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	ANKRD16 (ankyrin repeat domain 16) is a protein-coding gene.
Immunogen	Synthetic peptide conjugated to KLH, corresponding to C terminal amino acids 287-316 of Human ANKRD16(NP_001009943.1, NP_061919.1, NP_001009941.1)
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse, Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Size	100 µg
Buffer	Constituents: PBS
Preservative	0.09% Sodium Azide
Storage	Store at 4°C (up to 6 months). For long term storage store at -20°C

GENE INFORMATION

Gene Name [ANKRD16 ankyrin repeat domain 17 \[Homo sapiens \]](#)

Official Symbol	ANKRD16
Synonyms	ANKRD16; ankyrin repeat domain 16; ankyrin repeat domain-containing protein 16;
Entrez Gene ID	54522
Protein Refseq	NP_001009941.1
UniProt ID	Q6P6B7