



Rabbit Anti-Human CASP3 Polyclonal Antibody (DPABH-24873)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 55-160 of human CASP3 (NP_004337.2).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human, Mouse, Rat
Purification	Affinity purification
Conjugate	Unconjugated
Applications	WB, IHC, IF
Positive Control	Jurkat
Format	Liquid
Size	50 µl, 100 µl
Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Preservative	0.02% Sodium Azide
Storage	Store at -20°C. Avoid freeze / thaw cycles.

BACKGROUND

Introduction	This gene encodes a protein which is a member of the cysteine-aspartic acid protease
---------------------	--

(caspase) family. Sequential activation of caspases plays a central role in the execution-phase of cell apoptosis. Caspases exist as inactive proenzymes which undergo proteolytic processing at conserved aspartic residues to produce two subunits, large and small, that dimerize to form the active enzyme. This protein cleaves and activates caspases 6, 7 and 9, and the protein itself is processed by caspases 8, 9 and 10. It is the predominant caspase involved in the cleavage of amyloid-beta 4A precursor protein, which is associated with neuronal death in Alzheimer's disease. Alternative splicing of this gene results in two transcript variants that encode the same protein.

Keywords	CASP3; caspase 3, apoptosis-related cysteine peptidase; CPP32; SCA-1; CPP32B; caspase-3; CASP-3; CPP-32; apopain; procaspase3; protein Yama; PARP cleavage protease; cysteine protease CPP32; SREBP cleavage activity 1; caspase 3, apoptosis-related cysteine protease;
-----------------	--

GENE INFORMATION

Entrez Gene ID	836
-----------------------	---------------------

UniProt ID	P42574
-------------------	------------------------
