



Anti-HARBI1 (aa 11-40) polyclonal antibody (DPABH-24887)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Transposase-derived protein that may have nuclease activity (Potential). Does not have transposase activity.
Immunogen	KLH conjugated synthetic peptide between amino acids 11-40 from the N terminal region of Human C11ORF77 (NP_776172.1).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Size	100 µl
Buffer	Constituents: PBS
Preservative	0.09% Sodium Azide
Storage	Store at 4°C (up to 6 months). For long term storage store at -20°C

GENE INFORMATION

Gene Name	HARBI1 harbinger transposase derived 2 [Homo sapiens]
Official Symbol	HARBI1
Synonyms	HARBI1; harbinger transposase derived 1; C11orf77; putative nuclease HARBI1; harbinger transposase-derived nuclease;
Entrez Gene ID	283254
Protein Refseq	NP_776172.1
UniProt ID	Q96MB7
Function	metal ion binding; nuclease activity;