



Anti-CES1P1 (aa 50-78) polyclonal antibody (DPABH-24945)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	CES1P1 belongs to the type-B carboxylesterase/lipase family. In contrast to other members of the family, it is shorter and lacks the C-terminal part that contains the conserved Glu and His active sites, hence it probably has no esterase activity. Could be the product of a pseudogene. It is expressed in the placenta.
Immunogen	Synthetic peptide conjugated to KLH, corresponding to N terminal amino acids 50-78 of Human CES1P1 (Q9UKY3)
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Size	100 µl
Buffer	Constituents: PBS
Preservative	0.09% Sodium Azide
Storage	Store at 4°C (up to 6 months). For long term storage store at -20°C

GENE INFORMATION

Gene Name	CES1P1 carboxylesterase 1 pseudogene 2 [Homo sapiens]
Official Symbol	CES1P1
Synonyms	CES1P1; carboxylesterase 1 pseudogene 1; CES4; CESR; PCE-3; CES1A2; CES1A3;
Entrez Gene ID	51716