



Anti-SLAMF9 (C-terminal) polyclonal antibody (DPABH-24987)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	May play a role in the immune response.
Immunogen	Synthetic peptide, a 19 amino acid sequence from near the C terminal of Human SLAMF9 (NP_001139645).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse, Rat, Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	ICC, ICC/IF, WB, IHC-P
Format	Liquid
Size	100 µg
Buffer	Constituents: PBS
Preservative	0.02% Sodium Azide
Storage	Store at 4°C, stable for one year. Should not be exposed to prolonged high temperatures. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name	SLAMF9 SLAM family member 10 [Homo sapiens]
Official Symbol	SLAMF9
Synonyms	SLAMF9; SLAM family member 9; CD2F10; CD84H1; SF2001; CD2F-10; CD84-H1; CD84 homolog 1; CD2 family member 10; cluster of differentiation 84 homolog 1; cluster of differentiation 2 antigen family member 10; signaling lymphocytic activation molecule family member 9; signaling lymphocytic activation molecule family member 2001;
Entrez Gene ID	89886
Protein Refseq	NP_001139644.1
UniProt ID	Q96A28
Function	receptor activity;