



Anti-LRRC18 (aa 156-185) polyclonal antibody (DPABH-25221)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	May be involved in the regulation of spermatogenesis and sperm maturation.
Immunogen	Synthetic peptide corresponding to a region within C terminal amino acids 156-185 of Human LRRC18 (NP_001006940.3), conjugated to KLH.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse, Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Size	100 µg
Buffer	Constituents: PBS
Preservative	0.09% Sodium Azide
Storage	Store at 4°C (up to 6 months). For long term storage store at -20°C

GENE INFORMATION

Gene Name [LRRC18 leucine rich repeat containing 19 \[Homo sapiens \]](#)

Official Symbol	LRRC18
Synonyms	LRRC18; leucine rich repeat containing 18; UNQ933; UNQ9338; VKGE9338; leucine-rich repeat-containing protein 18;
Entrez Gene ID	474354
Protein Refseq	NP_001006940.3
UniProt ID	Q8N456