



Anti-CWF19L1 (aa 379-516) polyclonal antibody (DPABH-26114)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	CWF19L1 belongs to the CWF19 family.
Immunogen	Recombinant fragment, corresponding to a region within amino acids 379 and 516 of Human CWF19L1 (AAH08746).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Size	50 µl
Buffer	Constituents: 20% Glycerol, 0.1M Tris, 0.1M Glycine, pH 7.0
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name	CWF19L1 CWF19-like 2, cell cycle control (S. pombe) [Homo sapiens]
Official Symbol	CWF19L1
Synonyms	CWF19L1; CWF19-like 1, cell cycle control (S. pombe); CWF19-like protein 1;
Entrez Gene ID	55280
Protein Refseq	NP_060764.3
UniProt ID	Q69YN2
Function	catalytic activity; molecular_function;