



Anti-ENPEP (full length) polyclonal antibody (DPABH-26226)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Appears to have a role in the catabolic pathway of the renin-angiotensin system. Probably plays a role in regulating growth and differentiation of early B-lineage cells.
Immunogen	Recombinant full length protein corresponding to Human Aminopeptidase A.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, IHC-P
Format	Liquid
Size	50 µl
Buffer	Constituents: 10% Glycerol, 0.1M Tris glycine, pH 7.0
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name [ENPEP glutamyl aminopeptidase \(aminopeptidase A\) \[Homo sapiens \]](#)

Official Symbol	ENPEP
Synonyms	ENPEP; glutamyl aminopeptidase (aminopeptidase A); APA; CD249; gp160; glutamyl aminopeptidase; EAP; AP-A; aminopeptidase A; differentiation antigen gp160;
Entrez Gene ID	2028
Protein Refseq	NP_001968.3
UniProt ID	Q07075
Pathway	Metabolism of Angiotensinogen to Angiotensins; Peptide hormone metabolism; Renin-angiotensin system;
Function	aminopeptidase activity; metalloaminopeptidase activity; metalloaminopeptidase activity; metallopeptidase activity