



Anti-SYAP1 (aa 186-324) polyclonal antibody (DPABH-26437)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Immunogen	Recombinant protein fragment containing a sequence corresponding to a region within amino acids 186 and 324 of Human SYAP1 (AAH14657).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Size	50 µl
Buffer	Constituents: 10% Glycerol, 0.1M Tris, 0.1M Glycine, pH 7.0
Preservative	None
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.

GENE INFORMATION

Gene Name	SYAP1 synapse associated protein 2 [Homo sapiens]
Official Symbol	SYAP1

Synonyms	SYAP1; synapse associated protein 1; PRO3113; synapse-associated protein 1; SAP47 homolog; synapse associated protein 1, SAP47 homolog;
-----------------	---

Entrez Gene ID	94056
-----------------------	-----------------------

Protein Refseq	NP_116185.2
-----------------------	-----------------------------

UniProt ID	Q96A49
-------------------	------------------------
