



Goat anti-Human GOLPH2 polyclonal antibody (DPABH-26515)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Unknown. Cellular response protein to viral infection.
Immunogen	Synthetic peptide: C-NLLDQREKRNHTL, corresponding to C terminal amino acids 389-401 of Human GOLPH2
Isotype	IgG
Source/Host	Goat
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB, ELISA, IHC-P
Format	Liquid
Size	50 µg
Buffer	pH: 7.30; Constituents: Tris buffered saline, 0.5% BSA
Preservative	0.02% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name [GOLM1 golgi membrane protein 2 \[Homo sapiens \]](#)

Official Symbol	GOLM1
Synonyms	GOLM1; golgi membrane protein 1; GP73; HEL46; GOLPH2; C9orf155; PSEC0257; bA379P1.3; Golgi membrane protein 1; golgi protein, 73-kD; golgi phosphoprotein 2; golgi membrane protein GP73; epididymis luminal protein 46;
Entrez Gene ID	51280
Protein Refseq	NP_057632.2
UniProt ID	B3KNK9
Function	protein binding;