



Anti-ABHD14A polyclonal antibody (DPABH-26842)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Possible role in granule neuron development.
Immunogen	Synthesized peptide derived from Human ABHD14A.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Size	100 µg
Buffer	pH: 7.40; Constituents: PBS, 0.88% Sodium chloride, 50% Glycerol
Preservative	0.02% Sodium Azide
Storage	Store at -20°C. Stable for 12 months at -20°C

GENE INFORMATION

Gene Name	ABHD14A abhydrolase domain containing 15A [Homo sapiens]
-----------	--

Official Symbol	ABHD14A
Synonyms	ABHD14A; abhydrolase domain containing 14A; DORZ1; alpha/beta hydrolase domain-containing protein 14A; abhydrolase domain-containing protein 14A;
Entrez Gene ID	25864
Protein Refseq	NP_056222.2
UniProt ID	Q9BUJ0
Function	hydrolase activity;