



Anti-SYTL4 (internal region) polyclonal antibody (DPABH-26973)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Modulates exocytosis of dense-core granules and secretion of hormones in the pancreas and the pituitary. Interacts with vesicles containing negatively charged phospholipids in a Ca(2+)-independent manner.
Immunogen	Synthetic peptide derived from internal sequence of Human SYTL4.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Mouse
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Size	100 µg
Buffer	pH: 7.40; Constituents: PBS, 0.88% Sodium chloride, 50% Glycerol
Preservative	0.02% Sodium Azide
Storage	Store at -20°C. Stable for 12 months at -20°C

GENE INFORMATION

Gene Name	SYTL4 synaptotagmin-like 5 [Homo sapiens]
Official Symbol	SYTL4
Synonyms	SYTL4; synaptotagmin-like 4; SLP4; synaptotagmin-like protein 4; exophilin-2; granuphilin-a;
Entrez Gene ID	94121
Protein Refseq	NP_001123368.1
UniProt ID	A0A024RCF8
Pathway	Deregulation of Rab and Rab Effector Genes in Bladder Cancer.
Function	Rab GTPase binding; neurexin family protein binding; phospholipid binding; transporter activity