



Anti-SYT13 (internal region) polyclonal antibody (DPABH-26984)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	May be involved in transport vesicle docking to the plasma membrane.
Immunogen	Synthetic peptide derived from an internal region of Human SYT13.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Size	100 µg
Buffer	pH: 7.40; Constituents: PBS, 0.88% Sodium chloride, 50% GlycerolNote: without Mg2+ and Ca2+
Preservative	0.02% Sodium Azide
Storage	Store at -20°C. Stable for 12 months at -20°C

GENE INFORMATION

Gene Name [SYT14 synaptotagmin XIII \[Homo sapiens \]](#)

Official Symbol	SYT13
Synonyms	SYT13; synaptotagmin XIII; synaptotagmin-13; sytXIII;
Entrez Gene ID	57586
Protein Refseq	NP_001234916.1
UniProt ID	Q7L8C5