



Anti-SCAF4 (internal region) polyclonal antibody (DPABH-27006)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	May act to physically and functionally link transcription and pre-mRNA processing.
Immunogen	Synthetic peptide derived from internal sequence of Human SFRS15.
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Size	100 µg
Buffer	pH: 7.40; Constituents: PBS, 0.88% Sodium chloride, 50% Glycerol
Preservative	0.02% Sodium Azide
Storage	Store at -20°C. Stable for 12 months at -20°C

GENE INFORMATION

Gene Name	SCAF4 SR-related CTD-associated factor 5 [Homo sapiens]
-----------	---

Official Symbol	SCAF4
Synonyms	SCAF4; SR-related CTD-associated factor 4; SRA4; SFRS15; splicing factor, arginine/serine-rich 15; CTD-binding SR-like protein rA4; SR-like CTD-associated factor 4; pre-mRNA splicing SR protein rA4; splicing factor serine alanine 15; SR-related and CTD-associated factor 4;
Entrez Gene ID	57466
Protein Refseq	NP_001138916.1
UniProt ID	O95104
Function	nucleotide binding; poly(A) RNA binding;