



Anti-FAM201B (aa 144-176) polyclonal antibody (DPABH-27126)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	FAM201B (family with sequence similarity 201, member B) is a pseudogene.
Immunogen	Synthetic peptide conjugated to KLH, corresponding to a region within C terminal amino acids 144-176 of Human LOC339742 (Q86VD6).
Isotype	IgG
Source/Host	Rabbit
Species Reactivity	Human
Purification	Immunogen affinity purified
Conjugate	Unconjugated
Applications	WB
Format	Liquid
Size	100 µl
Buffer	Constituent: PBS
Preservative	0.09% Sodium Azide
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

GENE INFORMATION

Gene Name [FAM201B family with sequence similarity 202, member B \[Homo sapiens \]](#)

Official Symbol	FAM201B
Synonyms	FAM201B; family with sequence similarity 201, member B;
Entrez Gene ID	339742