



# Anti-TBB8B (aa 327-355) polyclonal antibody (DPABH-27137)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Immunogen</b>	Synthetic peptide conjugated to KLH, corresponding to a region within C terminal amino acids 327-355 of Human TBB8B (A6NNZ2).
<b>Isotype</b>	IgG
<b>Source/Host</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Purification</b>	Immunogen affinity purified
<b>Conjugate</b>	Unconjugated
<b>Applications</b>	WB
<b>Format</b>	Liquid
<b>Size</b>	100 µl
<b>Buffer</b>	Constituent: PBS
<b>Preservative</b>	0.09% Sodium Azide
<b>Storage</b>	Store at 4°C short term (1-2 weeks). Aliquot and store at -20°C long term. Avoid repeated freeze / thaw cycles.

## BACKGROUND

<b>Introduction</b>	Function: Tubulin is the major constituent of microtubules. It binds two moles of GTP, one at an exchangeable site on the beta chain and one at a non-exchangeable site on the alpha-chain.
---------------------	---

Similarities: Belongs to the tubulin family.

---

**Keywords**

LOC100510573; tubulin beta-8 chain B-like;

---